


# TINCTURE PRESS

This durable tincture press is lightweight, easy to use, quick to clean, and very low maintenance. Its compact size and light weight make it very portable and ideal for use in the field as well as at home, school, or at your business.



## CARE & MAINTENANCE

- 1) Wash all parts thoroughly before initial use. Subsequently, most of the time you will only need to wash the piston and pot after each use.
- 2) If you need to wash the clamp, make sure to lubricate screw threads with food grade oil. Run screw all the way in and out to prevent rusting.

## INSTRUCTIONS FOR USE

- 1) Pour off bulk of fluid.
- 2) Load herbs into pot.
- 3) Spread and tamp herbs evenly in pot.
- 4) Set piston on surface of herbs.
- 5) Insert screw into piston sleeve.
- 6) Tighten screw until snug.
- 7) Rest the clamp horizontally across a mixing bowl or half gallon measuring cup as the press continues to drain. See photo.
- 8) Every few minutes, you will be able to tighten the screw a little more and extract more fluid. (Tightening the screw incrementally rather than all at one time will extend the life of the clamp and may also yield more fluid.)
- 9) The key to getting the most fluid out of each batch is to allow a bit more time progressively between each cycle of tightening (e.g., tighten at 5 minutes, then 10 minutes, then 15 minutes, etc.)
- 10) Continue to tighten screw and drain until all fluid is extracted. This process will take half an hour to several hours depending on your herb.



## 10 YEAR LIMITED WARRANTY

Myronaturals warranties pistons and pots against defects and breakage within the scope of the intended use. Breakage due to dropage, use with a cheater bar, or any method not listed in the instructions for use section of this literature are not covered under our warrantee.

We appreciate your business and welcome your feedback.

HAPPY PRESSING!

[www.myronaturals.com](http://www.myronaturals.com)

541.302.2628